


Feces Collection Kit


MSK-FCK001&MSK-FCK002 is a medical device that is used to collect, preserve and transport feces samples for testing and research purpose, which plays a crucial role in detecting various medical conditions such as gastrointestinal conditions, monitoring digestive health, and conducting research on gut microbiota. The kit includes a feces collection tube, spoon or screw stick, samples collected are accurately and safely preserved for downstream analysis.



Collection Instructions

1. Take out the device from the box;
2. Unscrew the cap to collect the feces;
3. Put the sample into tube and screw it;
4. Write down the relevant information.

How to Use the Feces Collection Kit?



Product Features

- Physical performance is stable, with good sealing;
- Temperature tolerance ranges from -80 to 60°C;
- Non-invasive, easy collection, convenient transport;
- Two types available: Spoon and Screw sampling stick.

Product Pictures

