

MSK-HCB001 Sampling Kits have been mainly used for collection, stabilization and transportation of the cervical exfoliated cells before HPV screening. Normally, a sampling kit contains 1 cervical sampling swab, 1 collection tube with HPV transport medium, and 1 collection bag. We also provide HPV self-sampling solutions, customization is welcome.



HPV Sampling Kit Contains

Disposable Sampling Swab(s)

CBT003 is designed to collect cervical exfoliated cell, and they have been used for gynecological clinical diagnostics, the "tower" tip is mainly used for HPV screening. MandeLab has developed the cervical brushes which exceed the industry standard according to the ergonomic design. The material of the tip is flexible and harmless, thus it would provide safe and comfortable sampling method. Cervical brushes collect more samples and ensure that sufficient samples can be released into the transport medium.

HPV Transport Medium

The HPV transport medium is mainly composed of HBSS, Preservatives, Antibiotics, Ethylene Diamine Tetraacetic Acid, Tritonx-100 and Methanol. It is used to preserve and transport the collected cells, and only for in vitro analysis and testing, not for therapeutic use.

How to Use the HPV Sampling Kit?



Storage and Disposal

The optimum storage temperature is 2°C to 30 °C. Do not incubate, or freeze prior to use, improper storage will result in a loss of testing. Medical waste must be disposed of in compliance with the local legislation. Please take the appropriate precautions for infected material if necessary.

Product Pictures

