

TCT Sampling Kits

TCT Sampling Kits have been used for collection, stabilization and transportation of the endocervical exfoliated cells before TCT testing. Nowadays, many countries do HPV screening and TCT testing together, but some countries do it separately. Normally, a TCT sampling kit contains 1 cervical brush, 1 collection tube with cell transport medium, 1 instruction card and 1 collection bag. Customization is welcome, we can also use swab to replace cervical brush.



TCT Sampling Kit Contains


Cervical Brush

MS-CBE002 is a cervical brush used for collecting cervical exfoliated cell, it has been designed for TCT specimen collection, and Gynecological clinical diagnostics. The cervical brush is ideal for collecting large amount of cells and rapid elution of the specimens that instantly releases the cells into the transport medium. We have CE MDR-0197, FDA, TGA and MHRA.

Cell Transport Medium

It is mainly composed of Glucose, Preservative, Tris, Ethanol Absolute, Purified Water. TCT transport medium is used to preserve and transport the collected cells, and only for in vitro diagnostics, not for therapeutic use. We also provide HPV transport medium which could be used for both HPV &TCT testing.

How to Use the TCT Sampling Kit?



Storage and Disposal

The optimum storage temperature is 2°C to 30 °C. Do not incubate, or freeze prior to use, improper storage will result in a loss of testing. Medical waste must be disposed of in compliance with the local legislation. Please take the appropriate precautions for infected material if necessary.

Product Pictures

